

# MiSeq Sample Submission Guide

## Next-Generation Sequencing Service

Contact: Rhea Cornely  
Phone: 02 4570 1784  
Email: [ngs@westernsydney.edu.au](mailto:ngs@westernsydney.edu.au)

WESTERN SYDNEY  
UNIVERSITY



Research

## Metagenomics Sequencing Requirements

Please read the following submission guidelines carefully and prepare your samples accordingly.

**Failure to follow submission guidelines may result in the rejection of submitted samples or might incur additional costs.**

### Prior to submission

- Custom primer projects which require case-specific optimisation **must be discussed with the NGS Facility staff** prior to sample submission.
- Samples must be submitted **one month prior** to the scheduled library preparation start date to enable sufficient time for sample QC.
- **Amplicon samples** must be submitted on **separate plates** from any **genomic DNA samples** and the nature of the amplicons must be discussed with the NGS Facility staff.
- Download and fill in
  - MiSeq\_Sample\_submission\_sheet\_v2
- Fill out online
  - MiSeq Sample Submission Form
- Upload (online) your completed
  - MiSeq\_Sample\_submission\_sheet\_v2
- Print out to send with your DNA/amplicon samples your completed
  - MiSeq Sample Submission Form and
  - MiSeq\_Sample\_submission\_sheet\_v2

**All fields in the sample submission form must be completed. Incomplete submissions will not be accepted.**

## Format of submitted samples

The preferred format for all sample submissions are Eppendorf 96-well PCR plates. For small numbers of samples ( $\leq 24$ ), 1.5ml tubes are acceptable but 8-tube strips are preferred.

- $\leq 24$  samples may be submitted in either,
  - 1.5 mL LoBind tubes (Eppendorf Cat No.022431021) or
  - 0.2ml PCR tube strips (Eppendorf Cat No. 951010022)
  - in a fully skirted 96-well plate (Eppendorf Cat No. 0030128648) with flat caps (Eppendorf Cat No. 0030124847) or storage foil (Eppendorf Cat No. 0030127889)
- $\geq 25$  samples must be submitted in Eppendorf fully skirted 96-well plates (available from the WSU NGS Facility on request) following the plate layout detailed below with no well gaps.

**Please note: If wells gaps occur, each gap will be charged as a standard library preparation.**

A fee for sample transfer will be charged for  $> 24$  samples not submitted in Eppendorf fully skirted 96-well plates.

## Labelling

- Samples in 1.5 ml tubes must be labelled “customer initials - sample number” (eg BH – 1 or BH 1)
- Samples in plates must be labelled on the front rim of each plate with “customer initials – plate number – date” (eg BH – P1 – 20190616 or BH – Plate 1 – 16/6/2019)

## Required 96-well plate format

The plate layout must be followed to conform with our epMotion automated pipetting systems.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Sample 1	Sample 9	Sample 17	Sample 25	Sample 33	Sample 41	Sample 49	Sample 57	Sample 65	Sample 73	Sample 81	Sample 89
B	Sample 2	Sample 10	Sample 18	Sample 26	Sample 34	Sample 42	Sample 50	Sample 58	Sample 66	Sample 74	Sample 82	Sample 90
C	Sample 3	Sample 11	Sample 19	Sample 27	Sample 35	Sample 43	Sample 51	Sample 59	Sample 67	Sample 75	Sample 83	Sample 91
D	Sample 4	Sample 12	Sample 20	Sample 28	Sample 36	Sample 44	Sample 52	Sample 60	Sample 68	Sample 76	Sample 84	Sample 92
E	Sample 5	Sample 13	Sample 21	Sample 29	Sample 37	Sample 45	Sample 53	Sample 61	Sample 69	Sample 77	Sample 85	Sample 93
F	Sample 6	Sample 14	Sample 22	Sample 30	Sample 38	Sample 46	Sample 54	Sample 62	Sample 70	Sample 78	Sample 86	Sample 94
G	Sample 7	Sample 15	Sample 23	Sample 31	Sample 39	Sample 47	Sample 55	Sample 63	Sample 71	Sample 79	Sample 87	Sample 95
H	Sample 8	Sample 16	Sample 24	Sample 32	Sample 40	Sample 48	Sample 56	Sample 64	Sample 72	Sample 80	Sample 88	Blank

Please note: any well gaps will incur an additional charge.

## Preparing samples for submission

### Elution buffer

Please elute your samples in **nuclease-free water**. Elution of DNA in buffers other than water may necessitate a purification step\*. If this is not possible, a tube with the elution buffer has to be included in the submission, otherwise the submission will be rejected.

### Genomic DNA Requirements for sequencing run

- Please ensure that samples are of similar concentrations to enable optimal normalisation.
- It is preferable if all submitted samples are around 30-50 ng/ul.
- Environmental and plant samples can be discoloured which interferes with quantification measurements and downstream assays. A purification step\* will be added prior to processing these samples.
- The total DNA amount **must not be below 450 ng** (eg. 30ul of 15ng/ul DNA).
- If the DNA concentration is **< 15ng/μl**, but the total amount of DNA is sufficient, a purification and concentration step\* **will be done** before QC if all samples in the same plate have the same volume.

**\* A concentration and/or purification step will incur an additional charge /sample and is applied to all submitted samples. Please enquire for the price.**

Total amount	Concentration	Minimum Volume	A260/280	A260/230	Buffer
≥450 -1500 ng	15ng/μl – 150ng/μl	>30μl	1.8-2.0 <sup>^</sup>	>1.8 <sup>~</sup>	Nuclease- free water

### Customer prepared amplicons per sequencing run

- The Illumina overhang adapter sequence must be added to locus-specific sequences for the region to be targeted (please contact us for further guidance at [ngs@westernsydney.edu.au](mailto:ngs@westernsydney.edu.au)).

Input material	Concentration*	Minimum Volume	Purity	Buffer
≥150ng	≥10ng/μl	30ul	Purified PCR free of primers/primer dimers	Nuclease- free water

\* Concentration must be assessed using fluorescence-based quantification method (Qubit or Picogreen).

<sup>^</sup> A low 260/280 ratio (below 1.8) indicates protein contamination of nucleic acids.

<sup>~</sup> A low 260/230 ratio (below 1.8) indicates organic or carbohydrate compound contamination of nucleic acids.

### Quality Values for Metagenomic Samples

Metagenomic samples can often have A260/230 ratios lower than 1.8. Whilst this is not unusual with these sample types, it may be indicative of inhibitors being present in the samples and may prevent successful amplification. It is our standard procedure to flag poor A260/230 values to customers in their sample QC report and confirm that you wish to proceed with library preparation regardless of the indicated low 'quality' values.

### Supporting data

- If submitting concentration values based upon Nanodrop results, please submit the RAW DATA only – 260/280 and 260/230 values only (and do not complete the concentration column in the sample submission form).
- Gel images demonstrating successful PCR amplification of the submitted samples must be supplied with the sample submission.

## Shipping of samples

Ensure samples are shipped in leak-proof containers (either tightly-sealing tubes or well-sealed 96 well plates. Cross-contamination of samples can occur during shipping if containers are incorrectly sealed).

Please be aware that samples can be delayed in shipping and pack samples appropriately.

Please send your samples by courier to the address below:

Attention to: Dr Caroline Janitz  
 Next-Generation Sequencing Facility  
 Western Sydney University – Richmond Campus  
 Building L9, R1.14-1.18  
 50 Bourke St  
 Richmond, NSW 2753

# Metagenomic Primer Sets

## Bacterial Primers (16S)

341 Forward and 805 Reverse (Bacterial 16S 341F/805R)

16S variable regions v3 and v4

Name Sequence

341F 5' CCTACGGGNGGCWGCAG

805R 5' GACTACHVGGGTATCTAATCC

515 Forward and 806 Reverse (Bacterial 16S 515F/806R)

16S variable region v4

Name Sequence

515F 5' GTGCCAGCMGCCGCGGTAA

806R 5' GGACTACHVGGGTWTCTAAT

27 Forward and 519 Reverse (Bacterial 16S 27F/519R)

16S variable regions v1, v2 and v3

Name Sequence PCR

27F 5' AGAGTTTGATCMTGGCTCAG

519R 5' GWATTACCGCGGCKGCTG

799 Forward and 1193 Reverse (Bacterial 16S 799F/1193R)

16S variable region v5, v6 and v7

Name Sequence

799F 5' AACMGGATTAGATACCKG

1193R 5' ACGTCATCCCCACCTTCC



**CONSERVED REGIONS:** unspecific applications

**VARIABLE REGIONS:** group or species-specific applications

## Fungal Primers (ITS)

fiTS7 Forward and ITS4 Reverse (Fungal ITS2 fiTS7/ITS4)

amplifying ITS2

Name Sequence

fiTS7 5' GTGARTCATCGAATCTTTG

ITS4 5' TCCTCCGCTTATTGATATGC

ITS1oo Forward and ITS4 Reverse (Fungal ITS1 ITS1f/ITS2)

amplifying ITS1 and 2

ITS1oo FWD – 5' GGAAGGATCATTACCACA


ITS4 – 5' TCCTCCGCTTATTGATATGC

ITS1f Forward and ITS2 Reverse (Fungal ITS1&2 ITS1oo/ITS4)

amplifying ITS1

ITS1f – 5' CTTGGTCATTTAGAGGAAGTAA

ITS2 – 5' GCTGCGTTCCTCATCGATGC



## Protist Primers (18S)

18S forward (Euk\_1391f forward primer) and 18S (EukBr reverse primer) reverse amplifying SSU rRNA gene V9 hypervariable region

Name	Sequence
18Sf (Euk_1391f)	5' GTACACACCGCCCGTC
18Sr (EukBr reverse)	5' TGATCCTTCTGCAGGTTACCTAC